

New as of: **07.2017**



# USBCam4 Tether

User Manual

English



## Table of contents

1 - General Data .....	<a href="#">3</a>
1.1 - About these installation instructions .....	<a href="#">3</a>
1.1.1 - Structure of the document .....	<a href="#">3</a>
1.1.1.1 - Identification of the warning levels .....	<a href="#">3</a>
1.1.1.2 - Formats and symbols used. ....	<a href="#">4</a>
1.1.2 - Conventions .....	<a href="#">5</a>
1.1.3 - Contents and target group .....	<a href="#">5</a>
1.2 - General Information .....	<a href="#">6</a>
2 - Features and Requirements .....	<a href="#">7</a>
3 - User Interface .....	<a href="#">8</a>
4 - Installation .....	<a href="#">9</a>
4.1 - Confirm Functionality .....	<a href="#">12</a>
4.2 - Set Up USBCam4 with Sidexis 4 .....	<a href="#">14</a>
5 - First Use .....	<a href="#">16</a>

# 1 General Data

## 1.1 About these Installation Instructions

### 1.1.1 Structure of the document

#### 1.1.1.1 Identification of the warning levels

To prevent loss of data and /or failure of the software, please observe the warning and safety information provided in this document. This information is highlighted as follows:

#### **DANGER**

An imminent danger that could result in loss of data or failure of software.

#### **WARNING**

A possible dangerous situation that could result in the loss of data or failure of software.

#### **CAUTION**

A possible situation that could result in failure of software.

#### **NOTICE**

A possible situation that could result in misconfiguration of software.

#### **IMPORTANT**

Application instructions and other important information.

**Tip:** Information for facilitating work.

### 1.1.1.2 Formats and symbols used

The formats and symbols used in this document have the following meaning:

✓ Prerequisite 1. First action step 2. Second action step or → Alternative Action ↪ Result ➔ Individual action step	Requests you to do something.
<b>see “Formats and symbols used [→ 7]”</b>	Identifies a reference to another text passage and specifies its page number.
● List	Identifies a list.
“Command / menu item”	Identifies commands, menu items, or quotations.

### 1.1.2 Conventions

Example	Meaning
Click	A single click down and releasing of the left mouse button.
Double-click	Clicking and releasing of the left mouse button twice in quick succession.
"Ctrl+N"	On the keyboard: Press the <b>Ctrl</b> and <b>N</b> keys simultaneously.
Drag & Drop	Pulling along and letting go. Click the left mouse button on an element, hold it and then release the mouse button on a possible target.
Activate / deactivate a check box	Select / deselecting a check box with a mouse click in order to activate / deactivate its underlying function.
Activate / deactivate an option button	Selecting / deselecting an option button with a mouse click in order to activate / deactivate its underlying function.

### 1.1.3 Contents and target group

These installation instructions are intended for performing the installation of USB-Cam4 Tether and the use of USBCam4 Tether.

## 1.2 General Information

### During installation

#### **IMPORTANT**

During installation of any included software, no Windows updates should be running or be executed in the background as this may lead to problems.

### Existing software

If a previous version of USBCam4 Tether exists on the PC, it should be removed before starting the following installation process.

### Requirements

#### **WARNING**

This software requires that USBCam4 driver are installed.

## 2 Features and Requirements

USBCam4 Tether offers the ability to use your USBCam4's capture button inside Dentsply Sirona's Sidexis 4 software.


- Freeze video from USBCam4 inside Sidexis 4.
- Capture image from USBCam4 inside Sidexis 4.
- Adjust freeze length giving more time to decide on capture.
- Automatic detection of USBCam4 device.

Component	Minimum for 2D Workstations	Minimum for 3D Workstations
Operating System	Windows 7 Pro SP1 (32 or 64bit) Windows 8.1 Pro (64bit) Windows 10 Pro (64bit)	Please see Site Survey for information.
RAM	≥4GB	
Hard Disk	≥5GB	
	Dedicated Graphics Memory or Intel Onboard	<a href="http://DentsplySironaSupport.com">http://DentsplySironaSupport.com</a>
Graphics Card	Capture window indicator.	
DirectX	Capture timeout adjustment.	

Further information can be found on both [Sirona.com](http://Sirona.com) and [DentsplySironaSupport.com](http://DentsplySironaSupport.com) including installation quick guides, documentation and manuals.

Microsoft .NET 4.5 is required for USBCam4 Tether to operate. It is included with the installer.

### 3 User Interface



Component	Object	Function
A:	Button click indicator.	Green if the capture button is pressed.
B:	Sidexis 4 indicator.	Green if Sidexis 4 is open.
C:	Camera ready indicator.	Green is USBCam4 is connected.
D:	Capture window indicator.	Green if capture window is open.
E:	Capture timeout adjustment.	Controls how long video is frozen.
F:	Start with Windows checkbox.	If checked, USBCam4 Tether will start on boot.

## 4 Installation

USBCam4 Tether should be installed on any workstation that is intended to capture images with USBCam4.

Upon successful installation, it is recommended to store the USBCam4 Tether.msi in a safe location for future use.

### IMPORTANT

The following components are required for proper functionality.


Sidexis 4 Video Module

- Found in the Sidexis 4 installer under "Additional Components"


USBCam4 Driver

- Found at [SchickBySirona.com](http://SchickBySirona.com)


- 1 Locate the installer file "USBCam4 Tether.msi".
- 2 Double Click "USBCam4 Tether.msi".



3 Click "Next".




4 Confirm the installation directory, then click "Next".




NOTE: This installation directory should not be changed.

5 Click "Next".



6 Click "Close".





👉 The installation of USB Cam4 Tether is now complete.

## 4.1 Confirm Functionality


After successful installation, confirm functionality of USBCam4 Tether before attempting to use with Sidexis 4.

- 1 Locate USBCam4 Tether in the start menu.  
Start > All Programs > Sirona > USBCam4 Tether
- 2 Click "USBCam4 Tether"




☞ The USBCam4 Tether user interface is shown.

3 Connect the USBCam4 device.




↳ After a few seconds, the “Camera Ready” indicator turns green.

4 Press the capture button on the USBCam4 device.



↳ The “Button Click Detected” indicator should immediately turn green.

5 Open Sidexis 4.




↳ The “Sidexis 4 Detected” indicator should turn green.


## 4.2 Set up USBCam4 with Sidexis 4

It is sometimes necessary to configure Sidexis 4 to use the USBCam4 for capture.

- 1 Open Sidexis 4.  
Start > All Programs > Sirona > Sidexis 4
- 2 Open "Settings".



- 3 On the left, select the "Cameras" sub menu.




### IMPORTANT

If the "Cameras" sub menu is not present, please install the Video Module from the Additional components menu in the Sidexis 4 installer.

- 4 Select "Other"


- 5 In the list that appears, select “Schick Technologies USBCam4”.



## IMPORTANT

If “Schick Technologies USBCam4” is not in this list, please install the USBCam4 Driver from SchickBySirona.com

- 6 Click the “Show Video” button.




- ↳ If the USBCam4 is powered on, you should now see live video. At this point, everything is functional. Press the escape key on your keyboard to close this window.


## 5 First Use

If all previous steps are successful, continue by confirming functionality within Sidexis 4.


- 1 Open Sidexis 4.  
Start > All Programs > Sirona > Sidexis 4
- 2 Register a patient using “Register and Exposure”



- 3 Open the capture window by clicking the small camera button at the bottom.



- 4 With the capture window shown, power the camera on.  
You should now see live video.
- 5 Position the camera in the area you intend to capture.  
When ready, press the capture button on your USBCam4 1 time.  
The live video will freeze.




As the live video is frozen, you have approximately 3 seconds to decide on keeping the image or not. If you want to save this image, press the capture button 1 time within 3 seconds. If you do not want to keep this image, wait until the allotted time is expired, and the live video will return.

## NOTICE

This 3 second window for deciding to keep the image is a default time. It can be changed in USBCam4 Tether's user interface.

- 6 Capture an image with the method below.  
With live video displayed, press the capture button 1 time.  
Within 3 seconds, press the capture button again.  
 An image will be saved in the background.

- 7 You can continue to capture subsequent images in the same way.  
(Press once to freeze, press once more to capture - OR wait 3 seconds for another attempt)



- 8 Close the capture window with the method below.  
With live video being displayed, press the capture button on your USBCam4 twice within 1 second.



- 9 The capture window is now closed, displaying the images saved during the capture session.