

System Requirements SIDEXIS XG, Version 2.63

Requirements for the installation/operation of a SIDEXIS workstation

Operating Systems	Windows XP Service Pack 3 32 bit Windows Vista SP2 32 bit Windows 7 Professional SP1 32 bit (64 bit not tested) Windows 7 Ultimate SP1 32 bit and 64 bit Windows 7 Ultimate 64 bit (tested BootCamp 3.0.4 (325) on iMAC 24", Intel Core 2 Duo, 4GB, MAC OS 10.6.7) Windows 8.1 Professional 32 bit and 64 bit (tested Patch Level Dec. 09, 2014) Windows 10 Professional 32 bit and 64 bit (tested Patch Level March, 2016) Windows Server 2008 R2 Windows Server 2012 R2
Processor	32-bit (x86), min. 1GHz
Main memory	Min. 1GB
Graphic settings	1.024x768 pixels (recommended), 16.7 million colors ("TrueColor")
Graphics card	Min. 128MB
Hard disk	At least 5GB free hard disk space for program and 2D data. When using GALILEOS, please take note of the corresponding additional hard disk space required.
USB port	For Version 1.1 and higher, required for USB components only
Browser	Internet Explorer, Version 6.0 SP1 or higher
PDF software	Acrobat Reader, Version >= 8.0 (included on CD). Required for the PDF test report function.

Requirements for SQL Server

Operating systems	Windows Server 2008 SP2
	Windows Server 2008 R2 SP1
	Windows Server 2012
	Windows Server 2012 R2
	Windows XP SP3
	Windows Vista SP2
	Windows 7 Professional SP1
	Windows 8.1 Professional
Windows 10 Professional	
Processor	min. 1GHz
Main memory	Min. 1GB
Hard disk	2GB of available hard disk space
Browser	Internet Explorer, Version 6.0 SP1 or higher

Important - Installation .Net 3.5.x! Windows 8.1: Internet connection required / Windows Server 2012 / 2012 R2: Server Manager > Role Services and Features > Features > Add.

Important - Not released: Domain Controller!