

Compatibility Overview IO



Intraoral Sensor	Interface	Interface Firmware	PC SW	Sidexis 4 V4.3	Sidexis 4 V4.4	Operating Systems	
Xios XG Supreme / Schick 33	USB2	v2.38	IO sensor plugin V2.2	✓	✗	Win 10 x64	
		v2.38	IOSS V3.2	✗	✓	Win 10 x64 / Win 11 x64	
	USB3	v1.4.107	IO sensor plugin V2.2	✓	✗	Win 10 x64	
		v1.5.140	IOSS V3.2	✗	✓	Win 10 x64 / Win 11 x64	
	WiFi	v3.00.192	IO sensor plugin V2.2	✓	✗	Win 10 x64	
		v3.00.192	IOSS V3.2	✗	✓	Win 10 x64 / Win 11 x64	
	End of Sale						
	Xios XG Select / Schick Elite	USB2	v2.38	IO sensor plugin V2.2	✓	✗	Win 10 x64
v2.38			IOSS V3.2	✗	✓	Win 10 x64 / Win 11 x64	
WiFi		v3.00.192	IO sensor plugin V2.2	✓	✗	Win 10 x64	
		v3.00.192	IOSS V3.2	✗	✓	Win 10 x64 / Win 11 x64	
Xios AE / Schick AE	USB3	v1.4.107	IO sensor plugin V2.2	✓	✗	Win 10 x64	
		v1.5.140	IOSS V3.2	✗	✓	Win 10 x64 / Win 11 x64	
Xios Scan (end of sales 2020 ROW/21 JP)	N/A	v2.1.11	Xios Scan SW V2.0	✓	✓	Win 10 x64 / Win 11 x64	
Xios Plus (end of support July 2022)	Ethernet Wallbox	V02.12.00	IO sensor plugin V1.2.0	✓	✓	Win 8.1 x64 / Win 8 x64 / Win 10 x64	
3rd Party Systems							
Dürr Vista Scan	N/A	N/A	Dürr Vista Image Bridge V1.26.0.139	✓	✓	Win 10 x64	

Compatibility Overview EO



Device	Device Firmware	PC-SW	Sidexis 4 V4.3.1	Sidexis 4 V4.4	Windows 10/x64	Windows 11/x64
Axeos	06.01.01	2.0.1	✓	✓	✓	✗
	06.01.03 (SP3)	2.0.2 (SP3)	✓	✓	✓	✓
Orthophos SL 2D/3D	05.12.07	1.2.3	✓	✗	✓	✗
	05.15.01 (R3 SP1)	3.1.0 (R3 SP1)	✗	✓	✓	✓
Orthophos S 2D/3D	05.14.01	1.4.1	✓	✗	✓	✗
	05.15.01 (R3 SP1)	3.1.0 (R3 SP1)	✗	✓	✓	✓
Orthophos E	02.58.05	1.2	✓	✓	✓	✓
Orthophos XG 2D	02.58.05	1.2	✓	✓	✓	✓*
Orthophos XG 3D	04.14.07	Workstation: 1.5 RCU SW: 2.5.1	✓	✓	✓	✓*
Galileos Comfort/ Compact (End of support 30.11.2022)	4.14.07	Workstation: 2.3 RCU SW: 2.5.1	✓	✓	✓	✗
Galileos Comfort Plus	4.14.07	Workstation: 2.3 RCU SW: 2.5.1	✓	✓	✓	✗

* When changing the operating system to Windows 11, the functionality and security of the Orthophos XG product (XG 3, XG 5, XG Plus, XG 3D ready, XG 3D) is not negatively affected. Please continue to observe the explicit recommendations on IT and cybersecurity measures in the instructions for use.

SW Compatibility



Software	Sidexis 4 V4.2	Sidexis 4 V4.3.x	Sidexis 4 V4.4
Galileos Implant	✓	✓	✗
SICAT Suite V2.0.40	✗	✓	✓
X-Ray Control Book (legacy)	✓	✓	✗
X-Ray Control Book 4	✗	✗	✓
Sidexis iX	Sidexis iX Server V4.0 or higher	Sidexis iX Server V4.0 or higher	Sidexis iX Server V5.1.x
Windows 10 Pro /x64	✓	✓	✓
Windows 11 Pro /x64	✗	✗	✓

SiroCam Compatibility



SiroCam	Sidexis 4 V4.2	Sidexis 4 V4.3.x	Sidexis 4 V4.4
F	✓	✓	✓
AF	✓	✓	✓
AF+	✓	✓	✓
3	✓	✓	✗
C	✓	✓	✗
Digital	✓	✓	✗
UAF Plus	✓	✓	✗