

# Product Compatibility List

The following table provides compatibility information for Schick intraoral, panoramic, and X-ray source products.

Schick Hardware	Product Image	Compatible with Windows 32-bit OS	Compatible with Windows 64-bit OS
Schick AE		Yes	Yes
Schick 33		Yes	Yes
Schick WiFi		Yes	Yes
Schick Elite		Yes	Yes
Remote HS (retired)		Yes	Yes <sup>1</sup>
iPan HS Remote (retired)		Yes	Yes <sup>2</sup>
PanX USB (retired)		Yes	Yes <sup>3</sup>
Pan Elite (retired)		Yes <sup>4</sup>	No
iPan Standard Remote (retired)		Yes <sup>5</sup>	No

<sup>1</sup> Not supported on Windows 10.

<sup>2</sup> Not supported on Windows 10.

<sup>3</sup> Compatible when used with iPan HS Remote. Not supported on Windows 10.

<sup>4</sup> Supported on Windows XP only.

<sup>5</sup> Requires an Intel or VIA USB controller.

Schick Hardware	Product Image	Compatible with Windows 32-bit OS	Compatible with Windows 64-bit OS
SDX (retired)		Yes (up to Windows 7, 32-bit) <sup>6</sup>	No
USB Cam4 (retired)		Yes	Yes
USB Cam2 (retired)		Yes	Yes <sup>6</sup>
USB Cam (retired)		Yes (up to Windows 7, 32-bit) <sup>6</sup>	No
CDR 2000 (retired)		Yes (up to Windows 7, 32-bit) <sup>6</sup>	No
CDR Wireless (retired)		Yes (up to Windows 7, 32-bit) <sup>6</sup>	No

Schick Software	Windows 10	Windows 8 Pro	Windows 7 <sup>7</sup>	Windows Vista	Windows XP <sup>8</sup>	Windows Server 2010 <sup>9</sup>	Windows Server 2008 <sup>10</sup>
CDR DICOM 5	Yes <sup>11</sup>	Yes	Yes	Service Pack 2	Service Pack 3	Yes	Yes
CDR DICOM 4.5	No	Contact Support	Yes	Service Pack 2	Service Pack 3	Yes	Yes

Last updated on November 21, 2019

<sup>6</sup> Compatible when used with x64 Beta driver. Not supported on Windows 8.1 and Windows 10.

<sup>7</sup> Includes Home Premium, Business, and Ultimate versions.

<sup>8</sup> Extended support for Windows XP ended in 2014.

<sup>9</sup> SQL Server Express not supported on a Domain Controller.

<sup>10</sup> SQL Server Express not supported on a Domain Controller.

<sup>11</sup> Compatible with SQL Server Express 2012 and SQL Server Express 2014.