


Product Compatibility List

The following table provides compatibility information for Schick intraoral, panoramic, and X-ray source products.

Schick Hardware	Product Image	Compatible with Windows 32-bit OS	Compatible with Windows 64-bit OS
Schick AE		Yes	Yes
Schick 33		Yes	Yes
Schick WiFi		Yes	Yes
Schick Elite		Yes	Yes
Remote HS (retired)		Yes	Yes ¹
iPan HS Remote (retired)		Yes	Yes ¹
PanX USB (retired)		Yes	Yes ²
Pan Elite (retired)		Yes ³	No
iPan Standard Remote (retired)		Yes ⁴	No

¹ Not supported on Windows 10.

² Compatible when used with iPan HS Remote. Not supported on Windows 10.

³ Supported on Windows XP only.

⁴ Requires an Intel or VIA USB controller.

Schick Hardware	Product Image	Compatible with Windows 32-bit OS	Compatible with Windows 64-bit OS
SDX (retired)		Yes (up to Windows 7, 32-bit) ⁶	No
USB Cam4 (retired)		Yes	Yes
USB Cam2 (retired)		Yes	Yes ⁵
USB Cam (retired)		Yes (up to Windows 7, 32-bit) ⁶	No
CDR 2000 (retired)		Yes (up to Windows 7, 32-bit) ⁶	No
CDR Wireless (retired)		Yes (up to Windows 7, 32-bit) ⁶	No

Schick Software	Windows 10	Windows 8 Pro	Windows 7 ⁶	Windows Vista	Windows XP ⁷	Windows Server 2010 ⁸	Windows Server 2008 ⁹
CDR DICOM 5	Yes ¹⁰	Yes	Yes	Service Pack 2	Service Pack 3	Yes	Yes
CDR DICOM 4.5	No	Contact Support	Yes	Service Pack 2	Service Pack 3	Yes	Yes

Last updated on November 25, 2019

⁵ Compatible when used with x64 Beta driver. Not supported on Windows 8.1 and Windows 10.

⁶ Includes Home Premium, Business, and Ultimate versions.

⁷ Extended support for Windows XP ended in 2014.

⁸ SQL Server Express not supported on a Domain Controller.

⁹ SQL Server Express not supported on a Domain Controller.

¹⁰ Compatible with SQL Server Express 2012 and SQL Server Express 2014.