

Information on Software Version v2.0.37

Changes compared to v2.0.35

- **Plug and Play of USB cameras possible while Sidexis is running**

A USB camera - such as the SiroCam UAF Plus - can be plugged in with Sidexis already started and used directly if it was the last camera configured.

- **Product improvements**

This version of the video plugins fixes several bugs.

Compatibility and version history

Version	Contents / Update motivation (Changes compared to previous version)	Doc.- Chg. No.	TENEO software
v2.0.37	<ul style="list-style-type: none"> • Plug and Play of USB cameras possible while Sidexis is running • Product improvements 	130237	TENEO v8.x, 7.x, SINIUS v2.x, v1.x INTEGO v2.x, v1.x
v2.0.35	<ul style="list-style-type: none"> • Compatible with Sidexis 4 v4.3 and v4.2 		TENEO v8.x, 7.x, SINIUS v2.x, v1.x INTEGO v2.x, v1.x
v2.0.32	<ul style="list-style-type: none"> • Compatible with Sidexis 4 v4.2 • SIVISION 3 is no longer supported • New design • Saving still images without a patient being registered 	123815	TENEO v8.x, 7.x, SINIUS v2.x, v1.x INTEGO v2.x, v1.x

Technical documentation

Video Plugin Operators Manual, German (64 92 610 D3509)
 Video Plugin Operators Manual, English (64 93 204 D3509)
 Video Plugin Operators Manual, French (65 57 420 D3509)
 Video Plugin Operators Manual, Spanish (65 57 438 D3509)
 Video Plugin Operators Manual, Italian (65 57 446 D3509)